

# ISO 17025 Certified Calibration Tiles



ISO 17025 certified calibration tiles are essential tools for calibration laboratories and quality control facilities, providing the critical foundation for accurate gloss measurements across industries.

As the UK's premier manufacturer of glossmeters, Rhopoint Instruments offers ISO 17025 certified calibration tiles suitable not only for Rhopoint IQ and Novo-Gloss instruments but other gloss measuring instruments.

## Calibration companies trust our tiles to:

- Maintain the highest accuracy in their calibration services
- Provide documented traceability for quality audits
- Ensure consistent measurements across different facilities
- Meet international quality management requirements



## 100x100mm high gloss calibration standard with ISO 17025 certificate

| Mid gloss     | Code          |
|---------------|---------------|
| Mid gloss 20° | B1000-057/ISO |
| Mid gloss 60° | B1000-051/ISO |
| Mid gloss 85° | B1000-056/ISO |

| Low gloss     | Code          |
|---------------|---------------|
| Low gloss 20° | B1000-055/ISO |
| Low gloss 60° | B1000-052/ISO |
| Low gloss 85° | B1000-058/ISO |

100x100 tiles available at any GU value on request. Suitable for all gloss measuring instruments

## ISO 17025 calibration standards in aluminium holders for Rhopoint IQ and Novo-Gloss 2019 onwards

| Mid gloss     | Code          |
|---------------|---------------|
| Mid gloss 20° | B6000-055/ISO |
| Mid gloss 60° | B6000-051/ISO |
| Mid gloss 85° | B6000-056/ISO |

| Low gloss     | Code          |
|---------------|---------------|
| Low gloss 20° | B6000-057/ISO |
| Low gloss 60° | B6000-052/ISO |
| Low gloss 85° | B6000-058/ISO |

*IQ and Novo-Gloss tiles available.*



Rhopoint Instruments Ltd, St Leonards on Sea, UK • T: +44 (0)1424 739 622 • E: sales@rhopointinstruments.com • www.rhopointinstruments.com  
 Rhopoint Americas Inc., Troy, MI 48083, USA • T: 1.248.850.7171 • E: sales@rhopointamericas.com • www.rhopointamericas.com  
 Rhopoint Instruments GmbH, Weyarn, Germany • T: +49 8020 9214-988 • E: info@rhopointinstruments.de • www.rhopointinstruments.de